


Product Datasheet

Rat Neurokinin B ELISA Kit (orb566963)

Description	Rat Neurokinin B ELISA Kit
Reactivity	Rat
Range	78-5000pg/ml
Target	Tac3
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	P08435
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Rat
Tac3(Tachykinin-
3) ELISA Kit