


Product Datasheet

Rat Rrm2 ELISA Kit (orb566971)

Description	Rat Rrm2 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	Rrm2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q4KLN6
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Rrm2(Ribonucleoside-
diphosphate redu...