


## Product Datasheet

### Rat Cdk5r1 ELISA Kit (orb567047)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat Cdk5r1 ELISA Kit  |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 0.156-10ng/ml         |
| <b>Target</b>          | Cdk5r1                |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P61810</b>         |
| <b>Sensitivity</b>     | 0.094ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Rat  
Cdk5r1(Cyclin-  
dependent  
kinase 5 act...