


## Product Datasheet

### Rat Pfkp ELISA Kit (orb567070)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat Pfkp ELISA Kit    |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 0.312-20ng/ml         |
| <b>Target</b>          | Pfkp                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P47860</b>         |
| <b>Sensitivity</b>     | 0.188ng/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Rat Pfkp(6-phosphofructokinase type C) E...