


Product Datasheet

Rat Ghrl ELISA Kit (orb567114)

Description	Rat Ghrl ELISA Kit
Reactivity	Rat
Range	4.688-300pg/ml
Target	Ghrl
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q9QYH7
Sensitivity	2.813pg/ml
Expiration Date	Please enquire.



Rat
Ghrl(Appetite-
regulating
hormone)
EL...