


## Product Datasheet

### Rat Hsd3b1 ELISA Kit (orb567141)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat Hsd3b1 ELISA Kit  |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 78-5000pg/ml          |
| <b>Target</b>          | Hsd3b1                |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P22071</b>         |
| <b>Sensitivity</b>     | 46.9pg/ml             |
| <b>Expiration Date</b> | Please enquire.       |



Rat Hsd3b1(3  
beta-  
hydroxysteroid  
dehydro...