


Product Datasheet

Rat beta 1 Adrenergic Receptor ELISA Kit (orb567192)

Description	Rat beta 1 Adrenergic Receptor ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	ADRB1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P18090
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat
ADRB1(Adrenergic
Receptor Beta 1)
EL...