


Product Datasheet

Rat Alkaline Phosphatase, Intestinal ELISA Kit (orb567197)

Description	Rat Alkaline Phosphatase, Intestinal ELISA Kit
Reactivity	Rat
Range	0.781-50ng/ml
Target	Alkaline Phosphatase, Intestinal
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P15693
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Rat Alkaline Phosphatase, Intestinal(ALP...)