


Product Datasheet

Rat BMPR-II(Bone Morphogenetic Protein Receptor II) ELISA Kit (orb567237)

Description	Rat BMPR-II(Bone Morphogenetic Protein Receptor II) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	BMPR-II
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q91WY9
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat BMPR-II(Bone Morphogenetic Protein ...