


## Product Datasheet

### Rat CRT ELISA Kit (orb567308)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat CRT ELISA Kit     |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 62.5-4000pg/ml        |
| <b>Target</b>          | CRT                   |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P18418</b>         |
| <b>Sensitivity</b>     | 37.5pg/ml             |
| <b>Expiration Date</b> | Please enquire.       |



Rat  
CRT(Calreticulin)  
ELISA Kit