


Product Datasheet

Rat DPPIV(Dipeptidyl Peptidase IV) ELISA Kit (orb567343)

Description	Rat DPPIV(Dipeptidyl Peptidase IV) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	DPPIV
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P14740
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
DPPIV(Dipeptidyl
Peptidase IV)
ELI...