


## Product Datasheet

### Rat EGR1 ELISA Kit (orb567348)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat EGR1 ELISA Kit    |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 6.25-400pg/ml         |
| <b>Target</b>          | EGR1                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>P08154</b>         |
| <b>Sensitivity</b>     | 3.75pg/ml             |
| <b>Expiration Date</b> | Please enquire.       |



Rat  
EGR1(Early  
Growth  
Response  
Protein 1...