


## Product Datasheet

### Rat MIP-1beta ELISA Kit (orb567529)

|                        |                         |
|------------------------|-------------------------|
| <b>Description</b>     | Rat MIP-1beta ELISA Kit |
| <b>Reactivity</b>      | Rat                     |
| <b>Range</b>           | 15.625-1000pg/ml        |
| <b>Target</b>          | MIP-1 $\beta$           |
| <b>Storage</b>         | 2-8 °C for 6 months     |
| <b>Note</b>            | For research use only   |
| <b>Uniprot ID</b>      | <b>P50230</b>           |
| <b>Sensitivity</b>     | 9.375pg/ml              |
| <b>Expiration Date</b> | Please enquire.         |



Rat MIP-1 $\beta$ (Macrophage Inflammatory Prot...