


Product Datasheet

Rat PDGFsR-alpha ELISA Kit (orb567602)

Description	Rat PDGFsR-alpha(Platelet-Derived Growth Factor Soluble Receptor alpha) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	PDGFsR- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P20786
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat PDGFsR- α (Platelet-Derived Growth Fa...