


Product Datasheet

Rat PPAR-gamma ELISA Kit (orb567630)

Description	Rat PPAR-gamma(Peroxisome Proliferator-Activated Receptor gamma) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	PPAR-γ
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O88275
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat PPAR-
γ(Peroxisome
Proliferator-
Acti...