


## Product Datasheet

### Rat VASPIN ELISA Kit (orb567728)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rat VASPIN ELISA Kit  |
| <b>Reactivity</b>      | Rat                   |
| <b>Range</b>           | 93.75-6000pg/ml       |
| <b>Target</b>          | VASPIN                |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q8R4Z1</b>         |
| <b>Sensitivity</b>     | 56.25pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Rat  
VASPIN(Visceral  
Adipose  
Specific Ser...