


Product Datasheet

Rat alpha1-AGP ELISA Kit (orb567735)

Description	Rat alpha1-AGP(alpha1-Acid Glycoprotein) ELISA Kit
Reactivity	Rat
Range	6.25-400ng/ml
Target	α 1-AGP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P02764
Sensitivity	3.75ng/ml
Expiration Date	Please enquire.



Rat α 1-AGP(α 1-Acid Glycoprotein) ELISA...