


## Product Datasheet

### Rat beta-actin ELISA Kit (orb567739)

|                        |                          |
|------------------------|--------------------------|
| <b>Description</b>     | Rat beta-actin ELISA Kit |
| <b>Reactivity</b>      | Rat                      |
| <b>Range</b>           | 0.313-20ng/ml            |
| <b>Target</b>          | $\beta$ -actin           |
| <b>Storage</b>         | 2-8 °C for 6 months      |
| <b>Note</b>            | For research use only    |
| <b>Uniprot ID</b>      | <b>P60711</b>            |
| <b>Sensitivity</b>     | 0.188ng/ml               |
| <b>Expiration Date</b> | Please enquire.          |



Rat  $\beta$ -actin(Beta actin) ELISA Kit