


Product Datasheet

Rat ACTalpha2 ELISA Kit (orb567750)

Description	Rat ACTalpha2(Actin Alpha 2, Smooth Muscle) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	ACTα2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P62738
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
ACTα2(Actin
Alpha 2,
Smooth
Muscle)...