


## Product Datasheet

### Rabbit BDNF ELISA Kit (orb567926)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Rabbit BDNF ELISA Kit |
| <b>Reactivity</b>      | Rabbit                |
| <b>Range</b>           | 31.25-2000pg/ml       |
| <b>Target</b>          | BDNF                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Uniprot ID</b>      | <b>Q9BFJ8</b>         |
| <b>Sensitivity</b>     | 18.75pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Rabbit  
BDNF(Brain  
Derived  
Neurotrophic  
F...