


Product Datasheet

Rabbit IL-1alpha ELISA Kit (orb567970)

Description	Rabbit IL-1alpha(Interleukin 1 Alpha) ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/ml
Target	IL-1 α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P04822
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rabbit IL-1 α (Interleukin 1 Alpha) ELISA...