


## Product Datasheet

### Rabbit RUNX2 ELISA Kit (orb568010)

|                        |                        |
|------------------------|------------------------|
| <b>Description</b>     | Rabbit RUNX2 ELISA Kit |
| <b>Reactivity</b>      | Rabbit                 |
| <b>Range</b>           | 0.313-20ng/ml          |
| <b>Target</b>          | RUNX2                  |
| <b>Storage</b>         | 2-8 °C for 6 months    |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <b>Q6PKZ3</b>          |
| <b>Sensitivity</b>     | 0.188ng/ml             |
| <b>Expiration Date</b> | Please enquire.        |



Rabbit  
RUNX2(Runt-  
Related  
Transcription

...