


Product Datasheet

Rabbit TNF-alpha ELISA Kit (orb568030)

Description	Rabbit TNF-alpha(Tumor Necrosis Factor Alpha) ELISA Kit
Reactivity	Rabbit
Range	7.813-500pg/ml
Target	TNF- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P04924
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Rabbit TNF- α (Tumor
Necrosis
Factor
Alph...