


Product Datasheet

Rabbit alpha2-M ELISA Kit (orb568036)

Description	Rabbit alpha2-M(Alpha-2 Macroglobulin) ELISA Kit
Reactivity	Rabbit
Range	6.25-400ng/ml
Target	α 2-M
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	G1SQ70
Sensitivity	3.75ng/ml
Expiration Date	Please enquire.



Rabbit α 2-M(Alpha-2 Macroglobulin) ELIS...