


Product Datasheet

Rabbit LH ELISA Kit (orb568050)

Description	Rabbit LH ELISA Kit
Reactivity	Rabbit
Range	0.469-30ng/ml
Target	LH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Sandwich ELISA, Double Antibody
Uniprot ID	Q6IY74
Sensitivity	0.281ng/ml
Expiration Date	Please enquire.



Rabbit
LH(Luteinizing
Hormone)
ELISA Kit