


Product Datasheet

Sheep TNF-alpha ELISA Kit (orb568085)

Description	Sheep TNF-alpha(Tumor Necrosis Factor Alpha) ELISA Kit
Reactivity	Sheep
Range	15.625-1000pg/ml
Target	TNF- α
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P23383
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Sheep TNF- α (Tumor
Necrosis
Factor
Alpha...