


Product Datasheet

Sheep GHRL ELISA Kit (orb568093)

Description	Sheep GHRL ELISA Kit
Reactivity	Sheep
Range	31.25-2000pg/ml
Target	GHRL
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Uniprot ID	Q30DT1
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Sheep
GHRL(Ghrelin)
ELISA Kit