


## Product Datasheet

### Sheep GH ELISA Kit (orb568098)

|                          |                                    |
|--------------------------|------------------------------------|
| <b>Description</b>       | Sheep GH(Growth Hormone) ELISA Kit |
| <b>Reactivity</b>        | Sheep                              |
| <b>Range</b>             | 62.5-4000pg/ml                     |
| <b>Target</b>            | GH                                 |
| <b>Storage</b>           | 2-8 °C for 6 months                |
| <b>Note</b>              | For research use only              |
| <b>Application notes</b> | Sandwich ELISA, Double Antibody    |
| <b>Uniprot ID</b>        | <b>P67930</b>                      |
| <b>Sensitivity</b>       | 37.5pg/ml                          |
| <b>Expiration Date</b>   | Please enquire.                    |



Sheep  
GH(Growth  
Hormone)  
ELISA Kit