


## Product Datasheet

### Sheep ADH ELISA Kit (orb568110)

|                          |                                        |
|--------------------------|----------------------------------------|
| <b>Description</b>       | Sheep ADH ELISA Kit                    |
| <b>Reactivity</b>        | Sheep                                  |
| <b>Range</b>             | 15.625-1000pg/ml                       |
| <b>Target</b>            | ADH                                    |
| <b>Storage</b>           | 2-8 °C for 6 months                    |
| <b>Note</b>              | For research use only                  |
| <b>Application notes</b> | Competitive ELISA, Coated with Antigen |
| <b>Sensitivity</b>       | 9.375pg/ml                             |
| <b>Expiration Date</b>   | Please enquire.                        |



Sheep  
ADH(Antidiuretic  
Hormone) ELISA  
Ki...