

Product Datasheet

ALK antibody (orb5684)

Description	ALK antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ALK around the phosphorylation site of Tyr1278/1282/1283 (DI(p-Y)RAS(p-Y)(p-Y)RK)
Target	phospho-ALK (Tyr1278 + Tyr1282 + Tyr1283)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	174 kDa
Uniprot ID	Q9UM73
Dilution Range	Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.