


Product Datasheet

ALK (phospho-Tyr1604) antibody (orb5685)

Description	ALK (phospho-Tyr1604) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ALK around the phosphorylation site of Tyr1604 (GH(p-Y)ED)
Target	phospho-ALK (Tyr1604)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	174 kDa
Uniprot ID	Q9UM73
Dilution Range	Flow-Cyt=1µg /test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Jurkat
cells...