

Product Datasheet

Human DLL3 Protein (orb570269)

Description	Delta-like protein 3 (DLL3) is a transmembrane protein that belongs to the Delta/Serrate/Lag-2 (DSL)...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	DLL3
Preservatives	PBS with Arginine, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 50.4 kDa. The protein migrates as 50-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_058637.1
Purity	95%
Source	Human DLL3, His Tag (orb570269) is expressed from human 293 cells (HEK293). It contains AA Ala 27 - Leu 492 (Accession # Q9NYJ7-1).
MW	50.4 kDa
NCBI	NP_058637.1
Expiration Date	6 months from date of receipt.