


Product Datasheet

Human IL-31 RA Protein (orb570280)

Description	Interleukin-31 receptor subunit alpha is a protein that in humans is encoded by the IL31RA gene, als...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-31 RA
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 83.5 kDa. The protein migrates as 105-115 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-31 RA, Fc Tag (orb570280) is expressed from human 293 cells (HEK293). It contains AA Ala 20 - Glu 519 (Accession # Q8NI17-1).
MW	83.5 kDa
NCBI	XP_011541444.1
Expiration Date	6 months from date of receipt.



Human IL-31, His Tag capture on NTA-Seri...