


Product Datasheet

Human IL-23 R Protein (orb570284)

Description	Interleukin 23 receptor (IL-23R) is a type I cytokine receptor, and IL-23R pairs with the receptor m...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-23R
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 64.5 kDa. The protein migrates as 90-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human IL-23 R, Fc Tag (orb570284) is expressed from human 293 cells (HEK293). It contains AA Gly 24 - Gly 355 (Accession # Q5VWK5-1).
MW	64.5 kDa
NCBI	NP_653302.2
Expiration Date	6 months from date of receipt.



Captured
Human IL-23
R, Fc Tag
(Cat. No....)