


Product Datasheet

Human IL-23 R Protein (orb570286)

Description	Interleukin 23 receptor (IL-23R) is a type I cytokine receptor, and IL-23R pairs with the receptor m...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-23R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.9 kDa. The protein migrates as 55-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IL-23 R, His Tag (orb570286) is expressed from human 293 cells (HEK293). It contains AA Gly 24 - Gly 355 (Accession # Q5VWK5-1).
MW	39.9 kDa
NCBI	NP_653302.2
Expiration Date	6 months from date of receipt.



Biotinylated
Human IL-
23A & IL-12B
Heter...