

Product Datasheet

Human LRP-5 Protein (orb570288)

Description	Low-density lipoprotein receptor-related protein 5(LRP-5) is also known as BMND1, EVR1, EVR4, HBM, L...
Reactivity	Human
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	LRP-5
Preservatives	51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mlgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 96.6 kDa. The protein migrates as 95-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002326.2
Purity	95%
Source	Human LRP-5, Mouse IgG2a Fc Tag (orb570288) is expressed from human 293 cells (HEK293). It contains AA Glu 644 - Gln 1263 (Accession # O75197-1).
MW	96.6 kDa
NCBI	NP_002326.2
Expiration Date	6 months from date of receipt.