

Product Datasheet

Human PDGF R alpha / PDGFRA Protein (orb570289)

Description	PDGFRA(platelet-derived growth factor receptor A), also termed PDGFR α (platelet-derived growth facto...
Reactivity	Human
Endotoxins	1.0 EU per μ g
Conjugation	Unconjugated
Target	PDGF R alpha
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10 $\times$ His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 58.4 kDa. The protein migrates as 66-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_006197.1
Purity	95%
Source	Human PDGF R alpha, His Tag (orb570289) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Ala 528 (Accession # P16234-1).
MW	58.4 kDa
NCBI	NP_006197.1
Expiration Date	6 months from date of receipt.