

Product Datasheet

Rabbit PD-L1 / B7-H1 Protein (orb570290)

Description	Programmed cell death 1 ligand 1 (PDL1) is also known as B7-H, B7H1, MGC142294, MGC142296, PD-L1, PD...
Reactivity	Rabbit
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-L1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.1 kDa. The protein migrates as 32-44 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_008253343.1
Purity	90%
Source	Rabbit PD-L1, His Tag (orb570290) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - His 238 (Accession # G1SUI3-1).
MW	27.1 kDa
NCBI	XP_008253343.1
Expiration Date	6 months from date of receipt.