

Product Datasheet

Cynomolgus PVRIG Protein (orb570295)

Description	Cynomolgus PVRIG, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr 41 - Asp 171...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PVRIG
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.1 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_005549281.1
Purity	95%
Source	Cynomolgus PVRIG, Fc Tag (orb570295) is expressed from human 293 cells (HEK293). It contains AA Thr 41 - Asp 171 (Accession # A0A2K5WVV8-1).
MW	40.1 kDa
Uniprot ID	A0A2K5WVV8
NCBI	XP_005549281.1
Expiration Date	6 months from date of receipt.