

Product Datasheet

Mouse PVRIG Protein (orb570297)

Description	Human PVRIG (poliovirus receptor related immunoglobulin domain-containing protein), also known as CD...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PVRIG
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.6 kDa. The protein migrates as 48-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Ser 35 - Asp 165
Purity	95%
Source	Mouse PVRIG, Fc Tag (orb570297) is expressed from human 293 cells (HEK293). It contains AA Ser 35 - Asp 165 (Accession # A0A1B0GS01-1).
MW	40.6 kDa
Expiration Date	6 months from date of receipt.