

Product Datasheet

Human Thrombospondin-2 Protein (orb570301)

Description Thrombospondin-2 is also known as THBS2, TSP2, TSP-2. Thrombospondin-2 (TSP2) is an anti-angiogenic ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target Thrombospondin-2

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 129.7 kDa. The protein migrates as 180-200 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_003238.2


Purity 90%

Source Human Thrombospondin-2, His Tag (orb570301) is expressed from human 293 cells (HEK293). It contains AA Gly 19 - Ile 1172 (Accession # P35442-1).

MW 129.7 kDa

NCBI [NP_003238.2](#)

Expiration Date 6 months from date of receipt.



SDS-Page
analysis of
Human
Thrombospondi...