

Product Datasheet

SDHB Antibody (monoclonal, 2I3) (orb570315)

Description

SDHB Antibody (monoclonal, 213)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC, WB

Applications

Immunogen

E. coli-derived human SDHB recombinant protein (Position: A29-V280).

Form/Appearance

Lyophilized

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

29 kDa

Uniprot ID


P21912

Dilution Range


Western blot, 0.1-0.5µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of SDHB using anti-SDHB antibody.



IHC analysis of SDHB using anti-SDHB antibody.



IHC analysis of SDHB using anti-SDHB antibody.