

Product Datasheet

Metabotropic Glutamate Receptor 5/GRM5 Antibody (orb570385)

Description

Metabotropic Glutamate Receptor 5/GRM5

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, WB

Immunogen

E.coli-derived human Metabotropic Glutamate Receptor 5/GRM5 recombinant protein (Position: E24-D560).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

150 kDa

Uniprot ID


[P41594](#)

Dilution Range


Western blot, 0.5-1µg/ml, Mouse, Rat Direct ELISA, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.



Western blot analysis of GRM5 using anti...



Flow Cytometry analysis of U937 cells us...