



Product Datasheet

APG5L antibody (orb5704)


Description	APG5L antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human APG5L (101-200/275 aa)
Target	ATG5/APG5L
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	Q9H1Y0
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:100-500, Flow-Cyt=1ug/Test, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P analysis of human gastric using AT...



IHC-P analysis of Mouse brain using ATG5...



Western blot: Lane 1: Mouse Cerebrum; Lane

...