

Product Datasheet

RPL13A Antibody (orb570409)

Description

RPL13A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

E.coli-derived human RPL13A recombinant protein (Position:A2-V203).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

24 kDa


Uniprot ID
P40429
Dilution Range

Western blot, 0.1-0.25µg/ml, Human
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Direct ELISA, 0.1-0.5µg/ml, Human


Expiration Date

12 months from date of receipt.


1 2 3 4 5



WB analysis using anti-RPL13A antibody.L...



IHC analysis of RPL13A using anti-RPL13A...



IHC analysis of RPL13A using anti-RPL13A...