

Product Datasheet

BCL2 (Phospho-T56) antibody (orb570451)

Description	BCL2 (Phospho-T56) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Synthesized peptide derived from Human Bcl-2 around the phosphorylation site of T56.
Target	BCL2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10415
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:100-1:300, IP: 2-5 µg: mg lysate, IF: 1:200-1:1000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.