

Product Datasheet

JAK2 (Phospho-Y1007) antibody (orb570489)

Description	JAK2 (Phospho-Y1007) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthesized peptide derived from Human JAK2 around the phosphorylation site of Y1007.
Target	JAK2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60674
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:20000
Expiration Date	12 months from date of receipt.