

Product Datasheet

MCL1 (Phospho-S159) antibody (orb570727)

Description	MCL1 (Phospho-S159) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Synthesized peptide derived from Human Mcl-1 around the phosphorylation site of S159.
Target	MCL1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07820
Dilution Range	IF: 1:200-1:1000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.