

Product Datasheet

SGOL1 (Phospho-S14) antibody (orb570916)

Description	SGOL1 (Phospho-S14) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthesized peptide derived from Human Sgo1 around the phosphorylation site of S14.
Target	SGO1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5FBB7
Dilution Range	IHC-P: 1:100-1:300, ELISA: 1:5000
Expiration Date	12 months from date of receipt.