

Product Datasheet

MARK2 (Phospho-T596) antibody (orb570943)

Description	MARK2 (Phospho-T596) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Synthesized peptide derived from Human MARK2 around the phosphorylation site of T596.
Target	MARK2
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7KZ17
Dilution Range	IF: 1:200-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.